NMR-Spektroskopie Service

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NMR Spectroscopy Service

The Chemistry Department has **four nuclear magnetic resonance spectrometer (Bruker Avance 300, 360, 400 and 500 MHz)** on which liquid high-resolution as well as solid-state NMR spectra (*MAS NMR*) can be recorded.

To ensure an efficient NMR Service the following instructions should be followed:

- Generally, the NMR Service accepts only chemically pure substances. In order to estimate the required measuring time, the NMR Order form must be accompanied by a copy of a proton spectrum (including a parameter list) of the compound which should be examined.
 - Only orders that are accompanied by completely filled-in NMR Order forms will be processed by the NMR service.
 - For optimum results, the sample must be transferred to a 5 mm NMR tube of the quality UP (Ultra Precision). Furthermore, the sample must be filtered and well mixed to avoid concentration gradients. The height of the sample in the NMR tube must be at least 40 mm.
 - NMR measurements which are performed on the spectrometer **Bruker Avance III 500 MHz** must be approved by the signature of the Research Leader.
-